



PRESSALIT®  
CARE



R8468  
R8478  
RB2290

# Content

General safety . . . . .	4
Symbols used in this manual . . . . .	5
Specific safety . . . . .	5
Product label . . . . .	6
Intended use . . . . .	8
Intended operator profile . . . . .	8
Intended user profile . . . . .	8
Operating environment . . . . .	8
Assembly and installation . . . . .	9
Versions . . . . .	10
Operating instructions . . . . .	11
Height adjustment . . . . .	12
Cleaning . . . . .	13
Disinfection . . . . .	14
Maintenance . . . . .	15
Service inspection . . . . .	15
Trouble-shooting . . . . .	17
Technical data . . . . .	18
Disposal and recycling . . . . .	19

# General safety

## **WARNING**

- Read and understand the manual completely before use. Keep the manual close to the unit.
- Comply with this manual to avoid accidents involving serious personal injury.
- Never use or handle this product in other ways than are specified in this manual, as this can put your safety at risk and cause damage to the product.
- Make sure that everyone who is to install or use the product and its systems has the necessary information and access to this manual.

## **CAUTION**

- Please comply with these instructions to avoid the product becoming damaged or destroyed.
- Do not perform repairs, disassembly, or assembly operations, extensions, re-adjustments or modifications to this product. These must be carried out by Pressalit or by persons authorized by Pressalit only.
- Never allow children to operate this product unless under supervision. The product is not intended for play.
- Never use the product if it is defective or damaged. Do not use a defective product before it has been repaired.

## **NOTICE**

- The information in this manual is based on correct fitting of the product in accordance with our Assembly instruction.
- Pressalit Care cannot be held liable if the product is used in any way which differs from that described in the guide.
- We reserve the right to amend this manual without prior notice.

# Symbols used in this manual

The following symbols are used in this manual as well as the associated documentation supplied with the product.



**WARNING** indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.



**CAUTION** indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.



**IMPORTANT** indicates important information about handling and use of the product.



**NOTICE** is used to address practices not related to personal injury.



Pinching of hands. Take care to avoid injury to hands when in the vicinity of closing mechanical parts.



Electricity or potentially dangerous voltage.



Non-ionizing radiation.  
Electromagnetic compatibility warning.

## Specific safety



- Never leave a person lying on the bench without supervision.
- Always be careful when moving a person to and from the bench.






This product does not contain any of the hazardous substances mentioned in the Council Directive 2011/65/EC on RoHS.



Do not dispose as unsorted municipal waste. The product should be collected separately and returned to the designated recycling service.



Name and address of manufacturer.



Maximum user weight  
Safe working load  
Product weight

DD-MM-YYYY

Production date

XXXV

Mains voltage

IPXX

IP class

RZZZZYYYYXXX(X)

Article number

XX-XXHz


Frequency

XXW

Maximum power consumption

en

Product label placement ①.



## Intended use

This product is intended for bathing a person who is unable to stand up or sit down. The person is lying horizontally on the bench.

## Intended operator profile

The product is to be used by personnel trained in health care and trained in washing a person in their care.

## Intended user profile

The product can be used by all - adults, children, young, old, whose weight does not exceed the **maximum load of 175 kg**.

## Operating environment

For indoor use in a wet environment, such as a bathroom or wash room.

en

The product is intended to be used at room temperature, in a stand-alone environment, in private homes, nursing homes and hospitals under normal working conditions.

### NOTICE

If the product is installed in an atmosphere heavy in salt or chloride, the warranty of the nursing bench will be reduced from the normal 3 years to 1 year.

Salt and chloride have a corrosive effect on the products components and will cause a decrease in the life expectancy of the product. The damaging effects can, however, be reduced if the product is showered with tap water on a daily basis.

# Assembly and installation


For assembly and installation, please refer to the assembly instruction supplied with the product.

RB468 - RB478

**PRESSALIT**  
GROUP

---

XXX



- 📖 Assembly Instruction
- 📖 Monteringsvejledning
- 📖 Montageanleitung
- 📖 Notice de montage

- 📖 Montage-voorschrift
- 📖 Monteringsanvisning
- 📖 Monteringsveledning
- 📖 Instrucciones de montaje

Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.  
 Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.  
 Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.  
 Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.





Pressalit Oslo  
 Dr. Hovslor 17  
 0404 Oslo, NO

Telefon: +47 22 88 41 11  
 E-post: [oslo@pressalit.com](mailto:oslo@pressalit.com)  
 www.pressalit.com

RB2290

**PRESSALIT**  
GROUP

---

XXX



- 📖 Assembly Instruction
- 📖 Monteringsvejledning
- 📖 Montageanleitung
- 📖 Notice de montage

- 📖 Montage-voorschrift
- 📖 Monteringsanvisning
- 📖 Monteringsveledning
- 📖 Instrucciones de montaje

Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.  
 Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.  
 Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.  
 Selveien til utvalgte artikler finnes også på våre nettsider. Besøk vår nettside for mer informasjon.





Pressalit Oslo  
 Dr. Hovslor 17  
 0404 Oslo, NO


Telefon: +47 22 88 41 11  
 E-post: [oslo@pressalit.com](mailto:oslo@pressalit.com)  
 www.pressalit.com

## Versions

The nursing bench comes in two versions, and the safety rail in one version.



**R8468**



**R8478**



**RB2290**

# Operating instructions




Support the bench using your hand or body when folding it down, to avoid the risk of damage.

The built-in gas cylinders make it easy to fold the bench up against the wall after use.



**Take care to avoid injury to hands when in the vicinity of closing mechanical parts ①.**


**To pull the (optional) safety rail up:** Release the lock ② and pull the safety rail ③ from the underside of the bench, and swing it upwards. The safety rail will automatically lock into position.



**To fold the safety rail away:** Release the lock ② and swing the safety rail down to the underside of the bench. The safety rail will automatically lock into position.

en

# Height adjustment



When the bench is installed at the correct height (recommended installation height is 868 mm above the floor), it can be adjusted between 300 and 1000 mm above the floor.

The height adjustment is controlled by the wired remote control unit.

Press and hold the Up button ① to raise the bench, and press and hold the Down button ② to lower the bench.




**Make sure that there is no objects under the nursing bench when operating the height adjustment system.**

## Cleaning

Clean the nursing bench after each use. Use ordinary soapy water. Rinse off with fresh water and dry with a soft cloth.

The frame with the PVC mesh cover ① can be removed for cleaning.



The mesh cover can be zipped off and washed in regular washing machine at max 70° C.

Dry the water collection tray ②.

en

### IMPORTANT

- Do not use an autoclave or other type of steam cleaner.
- Avoid aggressive cleaning agents containing acids or chloride.
- Never spray the nursing bench directly with a high pressure cleaner.
- Avoid use of aggressive, abrasive, corrosive or chloride-based cleaning or disinfecting products.

# Disinfection

When using the nursing bench for multiple patients, it is necessary to both clean and disinfect the nursing bench between uses on different patients to prevent the risk of cross infection.

A list of suitable cleaning agents and disinfectants can be found at [www.pressaliticare.com](http://www.pressaliticare.com)



Always read the material safety data sheet (MSDS) and supplied instructions of the cleaning, disinfecting or de-scaling agent and comply with instructions for safe use, dilution, rinsing, drying, handling and storage of the agent.


# Maintenance

Expected service life of the nursing bench is 10 years, provided that maintenance is carried out as described below:

It is recommended on a weekly basis to run the motors on the product to their maximum positions by raising the product to its full height, and then lowering it to its lowest position, so that the internal parts are lubricated.

It is recommended that service inspection is done once a year.


## Service inspection




- Check the hinges ① on the bench and the safety rail ② . Tighten if necessary.
- Check the mesh cover ③ for damage and replace if necessary.



- Check the drain pipe connections ④, and check and clean the water trap (if fitted).




Remove the upper and lower retaining screws ⑤ on the lift unit cover and carefully remove the cover.



- Check the wall brackets ⑥ and tighten if necessary.
- Oil moving parts with WD40, especially the sliding brackets ⑦.
- Check all electrical connections to the control box ⑧.
- Check the chain wheel ⑨ for wear.
- Oil the chain ⑩.
- Put the unit cover back again and tighten the retaining screws.

# Trouble-shooting

Symptom	Cause	Action
The nursing bench does not stay in upright position.	One or more gas cylinders are defective.	Contact Pressalit A/S or your local dealer.
The product will not work, even though the remote control is used correctly.	Fault in power supply.	Check that all wires and plugs are intact and correctly connected. Contact Pressalit A/S or your local dealer.
The water does not run out of the water collection tray.	Adjustment screws too high or too low.	Adjust the adjustment screws ①.
The nursing bench runs up and down unevenly.	The sliding brackets are not properly lubricated.	Oil the sliding brackets as described under Maintenance,



Contact Pressalit Care or your dealer if the problem cannot be remedied. Do not attempt to repair the product yourself, as this will invalidate the guarantee.

## Technical data

Max load	175kg
Materials	Powder laquered stainless steel Plastics PVC (mesh cover)
Dimensions (width x length)	R8468: 875 mm x 1800 mm
	R8478: 899 mm x 1800 mm
Weight	R8468: 67 kg
	R8478: 73 kg
IP class	The product has an ingress protection class rating of IPX6. The product may not be immersed in water, but is designed to withstand powerful jets (12.5 mm nozzle) of water from any direction.
<b>Power specifications</b>	
Max consumption	70 W
Voltage	Mains connection: 230V/120V/100V
Frequency	50/60 Hz
Current	3.7 A

# Disposal and recycling

The unit contains reusable materials. There are no known hazards or risks associated with disposal of the product. All components can be discarded after having been cleaned and disinfected.

When disposing of the product, we recommend that it be disassembled and broken down into different waste groups for recycling or combustion.

The following table provides information on the recycling and handling of the product.

<b>Recyclable components</b>	<b>Material</b>
Steel	Steel can be recycled by being re-melted and included as a secondary material in the production of new steel.
PVC and polystyrene	Hard PVC is granulated and re-melted. Use plastic sorting and recycling systems where they are locally available, otherwise incinerating plastic has the added benefit of generating energy, which can be used for, for example, combined heat and power production.
Electronics	The electronic components should be collected separately and returned to the designated recycling service.