



*Access Your World*

## ***Linido / Handicare***

# **Folding Changing Table User Manual**

**(Models Li2403.1901-02;Li2590.1251-02 &  
Li2403.1501-02;Li2590.1251-02)**

### **Important safety information**

- Please read this manual thoroughly. We take no responsibility for any unspecified uses.
- The information in this manual pertains to correctly installed equipment.
- The product is not intended to be played with.
- Do not use the changing table if it is out-of-order or damaged.
- This manual is subject to alteration without notice.
- Support the table while unfolding to avoid damage.

### **Using the table**

The table can fold against the wall and be secured with the provided bracket. It can also hang from the mounting hardware in a downward position.

Unfold the two leg supports where they hinge to the bed to use the table for showering or changing. **Push both leg supports completely into the black hinges.** This locks the leg supports into position, providing safety and stability. Both legs should rest securely on the floor.

**Make sure both leg supports are fully disengaged** from the black hinges before folding the table to return it to the wall for storage. If the leg supports are not unlocked before being folded to the table, the hinge assembly may be damaged. Secure the folded legs against the table with the provided clips.

**Do not leave anyone unattended while they are occupying the table.**



## Cleaning

All components are water compatible. Most general purpose cleaners and disinfectants can be used. Avoid harsh brushes, abrasives, undiluted bleach products, or cleaning agents containing acids or chloride. Rinse off with clear water and dry with a soft cloth. Clean after every use.

## General Maintenance

We recommend preventative maintenance every three months:

- Check the table folding operation.
- Check all fasteners and anchors, and tighten as necessary.
- Check stretcher bed for any damage and clean as recommended.
- Report damaged or missing parts.
- No lubrication or adjustment is necessary.

Verify your purchase details before carrying out repairs so you don't void your warranty:

*We guarantee all product components against manufacturing defects for five years, and parts repair/labor for one year, from the invoice date. In most instances the affected parts are replaced, not repaired.*

*Parts damaged due to improper product assembly, installation, or procedure; lack of recommended maintenance including cleaning; or lost after delivery cannot be included in this warranty.*

For any questions or problems, please contact us at 1 800 577-1555 or [www.max-ability.com](http://www.max-ability.com).